

croquis



CROQUÍIS®

Samyang Biopharmaceuticals Corporation

Business Areas	4
Major Achievement	
Major Locations	5

SYB Aesthetics

Art Technology	6
----------------------	---

croquíis®

Brand Motive	8
MONO	10
SCREW	12
BARBED CUTTING	14
BARBED MOLDING	16

Guaranteed Quality

Stress Test	20
BHS Test	22

* BHS (Barb Holding Strength)

Brand Vision

croquíis® as Market Standard	24
------------------------------------	----




Samyang Biopharmaceuticals Corporation


Through accumulated experience and technology

With more than 20 years of experience in suture manufacturing and polymer technology, Samyang Biopharm's medical devices business intends to further expand its portfolio of wound management, devices for minimally invasive surgery and medical aesthetics products

Business Areas



Major Achievements



Major Location



Samyang Discovery Center
(Headquarter & R&D Center)

Completion	2016
Related Department	Marketing & Sales / Regulatory Affairs / Planning / R&D Center
Location	Seongnam-si, Gyeonggi-do, Korea

Pharmaceuticals Plant



COMPLETION
1997

PRODUCT
Injection, API

LOCATION
Daejeon, Korea

CERTIFICATIONS
EU GMP, PMDA (Japan), GMP, KGMP (Korea)

Medical Devices Plant




COMPLETION
1996

PRODUCT
Sutures, Hemostats, Patches

LOCATION
Daejeon, Korea

CERTIFICATIONS
KGMP, CE Mark, ISO13485 (Medical Devices), KGMP (Pharmaceuticals)

Samyang Biopharm USA Inc.



COMPLETION
2018

BUSINESS MODEL

- NRDO (No Research Development Only)
- Therapeutics for Oncology and Rare Diseases
- License-in / Out


LOCATION
Cambridge, MA, USA



SYB Aesthetics

Everyone is born beautiful. The definition of beauty has changed from wanting to look like someone to bringing out your beauty within. Samyang believes such value as art. SYB Aesthetics is derived from combining art with technology. By bringing out your embedded beauty with Samyang's technology, SYB Aesthetics provides a helpful and safe solution to stay young and beautiful.


Art Technology





croquis®

Brand Motive

croquis® is motivated by the drawing technique of describing an object in a short period of time. Find your embedded beauty through skin rejuvenation and tissue repositioning.



RETURNING TO THAT YOUTH



MONO

Tightening ●●● | Volume ●○○

croquis MONO stimulates neocollagenesis, which not only improves your overall skin quality but also gets completely absorbed into the body over time.

Material : PDO (Polydioxanone)

Contains : 60ea / box (10ea / Al pouch * 6)

Thread Type : _____

Gauge	Cannula		Thread	
	Length (mm)	Length (mm)	Length (mm)	USP
25G	60	90	150	5 - 0
	90	90		
27G	60	90	50	6 - 0
	50	70		
29G	40	50	50	7 - 0
	25	30		
30G	40	50	30	7 - 0
	30	30		
31G	30	30	50	
	40	50		





Tightening ●●● | Volume ●●●

croquis[®] SCREW stimulates neocollagenesis, which improves your overall skin and wrinkle through volumization effect. Likewise, SCREW type also gets absorbed into the body over time.



Material : PDO (Polydioxanone)
 Contains : 60ea / box (10ea / Al pouch * 6)
 Thread Type :

Cannula		Thread	
Gauge	Length (mm)	Length (mm)	USP
26G	50	70	5 - 0
	60	90	
	90	130	
27G	40	52	6 - 0
	50	70	
	90	130	
29G	40	52	6 - 0
	50	70	

B RINGING THAT SMOOTH AND TIGHT SKIN



PULLING AND HOLDING THAT BEAUTY



BARBED CUTTING

Skin Quality ●●● | Lifting ●●●

croquis[®] BARBED CUTTING is used for tissue fixation and repositioning. With absorbable component and neocollagenesis, barbed cutting type provides a younger and beautiful look.

Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / Al pouch * 5)


Thread Type : 

product	Cannula		Thread		Barb
	Gauge	Length (mm)	Length (mm)	USP	
A Type Barb (Needle)	23G	90	150	4 - 0	 (Bi-directional)
		60	110	3 - 0	
B Type Barb (Needle)	23G	90	150	4 - 0	 (3D)
		60	110	3 - 0	
B Type Barb (L-Cannula)	20G	90	140	2 - 0	
		60	105	0	
		18G	100	150	2

 **BARBED MOLDING**


Skin Quality ●●● | Lifting ●●●●+

croquis[®] BARBED MOLDING is used for a stronger lift. By using the pine tree-like barb shape, barbed molding type gives a helpful solution to skin laxity.

Material : PDO (Polydioxanone)
 Contains : 30ea / box (6ea / Al pouch * 6)
 Thread Type : 

Cannula		Thread	
Gauge	Length (mm)	Length (mm)	USP
18G	100	185	2
19G	100	185	1

EN HANCING AND KEEPING THAT LIFT



Guaranteed Quality

High stability of croquis can solve quality-related issues during storage and distribution processes. Such stability provides mutual trust to global customers.

Storage Test




TYPE : MONO
USP : 6 / 0

Product	croquis MONO	Competitor A	Competitor B
Type	MONO	MONO	MONO
USP	6 / 0	6 / 0	6 / 0
Tensile Strength 5 years simulation (%)	75% Strength Remained	Not Testable (Powdered)	Not Testable (Decomposed)

Accelerated Aging Condition : 5 year-storage simulation by Q10 Regulation (Ref. ASTM F1980-16)

Result after 5-year accelerated aging test with competing brands




BHS Test



BHS : Barb Holding Strength (kgf)
USP : 0


croquis BARBED CUTTING type's high BHS demonstrates effective facial lifting.



Brand Vision

croquis® as Market Standard

SYB Aesthetics provides high and safe quality products through well-controlled and strict manufacturing process.



Category	Description
Pouch Characteristic	Easy Peel
Sterilization	EO Gas / Drying
Storage Period	3 Years
Storage Condition	Room Temperature

Draw your beauty,
Capture the moment

croquis®

